

Welcome to an Educational break

April 2020

THIS DOCUMENT IS PUBLIC

April 2020
Air Liquide Group



To start



Mute the microphone if you don't speak ... and remember to switch it on when you speak in Q&A!



Turn off the camera to save bandwidth



Consider using your headphones!



Adelino Fernandes

**South West Europe
Research & Analysis
Offer Deployer Manager**



Experience @AL: 12 years in Air Liquide: mainly in Specialty Gases Sales & Marketing (Project, Product and Market Manager) for Portugal and Spain. Technical Career Leader Senior Expert. Business Lab Pioneer and Business Developer for Isotopes & Natural Air WW.

Background: MSc. in Mechanical Engineer, BMus in music, Technical vehicle inspector, Research in High-Speed Machining and Hydrolift Maintenance, 4 patents register.

Contents

Presentation



20 min

**Air Liquide
in Brief**

**Gas cylinder
Handling**

Safety

**Building
the Future**



15 min

**Questions
& Answers**

1

Air Liquide in Brief

2019 Key Figures



~67,000
EMPLOYEES



PRESENT IN
80 COUNTRIES



MORE THAN
3.7 MILLION
CUSTOMERS & PATIENTS



REVENUE
€22bn



NET PROFIT
(GROUP SHARE)
€2.24bn



INVESTMENT
DECISIONS
€3.7bn



300 million
of innovation expenses/year

An Innovative Group

Innovation is at the heart of the Group's customer-centric transformation strategy



330
new patents
filed in 2019

4,300
employees⁽¹⁾
contribute to
innovation

€317m
innovation
expenses⁽¹⁾

(1) 2019 Figures OECD Definition.

THIS DOCUMENT IS PUBLIC

AIR LIQUIDE, A WORLD LEADER IN GASES, TECHNOLOGIES AND SERVICES FOR INDUSTRY AND HEALTH

2

Gas cylinder Handling

Cylinder main elements...



C

CAP

Used to move the cylinder and Protect the Valve.



V

Valve

The most fragile part in the cylinder. Always close after using.



B

Body

Cylinder that contain high pressure gas.

Gas color code



light Blue
Oxidizing



maroon
Acetilene



color mix
variety risk



white
Oxygen



deep green
Inert



blue
Protoxide



red
flammable



yellow
**toxic
corrosive**

and pictogram



- Gas under **pressure**: danger of explosion if heated;
- They are inert;
- Asphyxiant** at high concentrations.



- Oxidizing** gases: Can cause or aggravate a **fire**;
- Keep valves and fittings **free of oil and grease**.



- Flammable** gases: Can be easily ignited.
- Keep away** from heat, hot surfaces, sparks, flames and all other **sources of ignition**.

Cylinder handling PPE



P

The use of a **safety cap** or **helmet** will be recommended or required depending on the working environment and the specific risk assessment.

P

The use of **gloves** and **long sleeves** is recommended and for the use of oxidizing or flammable gases the clothing should be **antistatic and fireproof**.

E

The use of **safety shoes** are mandatory.

3

Safety

Disclaimer: It is a knowledge refreshment not replacing a formal training

Main Risks...

When handling a cylinder:

**Cylinder
Fall**

**High
Pressure**

**Gas
Composition**

**Technical
Unknown**

Due to low stability

Up to 7° of
inclination
without falling.

By Shocks

Not attached
to anchorage
point.

By dragging them...

Regulator

- Bad connection;
- Inappropriate
regulator;
- Wrong purge.

Materials compatibility

- Improper use;
- Regulator material
incompatibility;
- Handle high pressure.

Contempt or ignorance risk

- Operation by untrained persons;
- Improper Use.

Confined Spaces, something to take attention...

If you...

- hear a gas leak;
- have symptoms of general malaise;
- see a colleague is unconscious;
- see an oxygen detector with an alarm signal;
- work in areas with inert gases and/or gases that displace oxygen or if you work in areas where gases are used.

Applies Safety Procedure


Consider the risk of ASPHYXIA

**Think before
you act**



**Leave the
Danger Area**

Anoxia...



Purge...

Why is it necessary to purge the equipment?

Regulator

Pipeline

Equipment



Contamination

Oxygen, humidity, dust...

Chemical reactions




O₂ concentration
(ppm)

Purge of nitrogen tubing to remove
oxygen in air

(N₂: 0,2 l/min, 8 bar, 20°C)

1 ppm



4

Building the Future

Our innovation approach

LABTOP™ is an advanced built-in regulator with on/off handwheel:

**Work
Safely**

**Protecting &
Ergonomic
CAP**

**Easy
Operations**

**Gain
Performance**

**SAVE:
Gas, Time,
Money**

1



outlet raccord:
6 mm and 1/8"



2



3



4



5



Increase your lab efficiency

S

Working Safely

- Quick **shut-off** of the gas;
- User is **not exposed** to high-pressure;
- Risk of **leakage and intoxication** is **reduced** significantly.





Increase your lab efficiency

P

Gain Performance

- **Non Residual Valve:** prevents contamination guaranteeing an qualitative calibration;
- The **double-stage regulator** ensures stability of the outlet pressure needed in calibration.



Increase your lab efficiency

S

Save: Gas, Time and Money

- **Prevent expensive losses** and leaks (On/Off);
- Better manage/**reduction your stock** (gauge);
- **Quickly start** your operation (connector);
- **Save you regulator and maintenance costs.**

Where can I use my LABTOP cylinder

Mixtures

Gaseous with non corrosive molecules

LABTOP
user guide



Digital
CoA



Calibration

Mixtures



Pressure
Adjustable from
0 to 10 bar

Flow rate
1 to 10 L/min

**Double Stage
regulator**
incorporated.

**More stable
cylinder**

Thank you for being with us!

April 2020

**Open session for
Questions & Answers...**

More details can be found in your email box.

THIS DOCUMENT IS PUBLIC

April 2020
Air Liquide Group

OBRIGADO

MERCI

GRACIAS

GRAZIE

THANK
YOU

